

DISPOSABLE DRESSING SET BENEFITS



- ✓ Components are presented in a single, sealed, sterile pack.
- ✓ Storage needs, set up time and waste are all reduced by use of one pack instead of several items
- ✓ Procedures and working practices can be standardised throughout the trust.

Disposable Dressing Set

FOR MEDICAL PROFESSION AND HOME USE

Disposable Dressing Set is specialized in custom made sets and procedure trays. We offer you a wide range of components to choose.

Manufactured and distributed by

GHM
Good Health Medical Co.,Ltd.



GOOD HEALTH MEDICAL CO.,LTD.

11/25 Moo4, Kanjanapisek Road
Sao-Tonghin, Bangyai
Nontaburi 11140 Thailand

E-mail: info@ghm.co.th Fax: 02-967-8086

Sales representative: +66626985428,
+66818150719,+66993968856

PRODUCTS SPECIFICATION



Medical Gauze Swab

Gauze: 3x3 8ply, Yarn: 40's, Mesh 24x20 made from pure 100% cotton, high absorbent, softness

Wrapping drape

Instrument Wrapping drape 40x40cm made of natural fiber-based materials, each sheet is biodegradable and no emit odors, toxic waste or other residues in a incinerator



Fenestration Drape

Fenestration Drape 60x60cm, central rectangular 10.5x14cm with tape, hole diameter 10cm, Made from Blue color Polyesthylene (PE) and Polypropylene (PP) 55gsm



Plastic Cup

Plastic Cup 30ml made from medical grade Polyethylene terephthalate (PET)



Plastic Forceps

Plastic Tooth Forceps made from Polypropylene (PP) medical grade, size 5"



Plastic Tray

Plastic Tray with 3 compartments made from medical grade Polyethylene terephthalate (PET)



Absorbent Cotton Balls

Cotton Balls: 100% Pure Nature Cotton Clean White, 0.35g Soft, Smooth and Non Fluorescent



Surgical Glove

Latex Examination Glove Power Free, Food grade according to EN71 part 3 : 1994+ A1: 2000 + AC:2002




Plastic Tray

made from medical grade Polyethylene terephthalate (PET) Wide 27cm x High 3cm



Sponge handle

Plastic handle medical sponge, length 16cm, sponge size 5cm.



Blade #11

Disposable plastic handle blade No. 11

UNIVERSAL DRESSING SET

Hemodialysis Dressing Set DLC

- 1 x Tray of 3 Compartments
- 1 x Plastic forceps
- 1 x Drape blue wrap 40 x 40cm sterile field
- 1 x Drape blue napkin 40 x 40cm
- 8 x Gauze Swab 3 x 3 8ply
- 8 x Absorbent Cotton Balls
- 2 x Glove S



Hemodialysis Dressing Set AVF

- 1 x Tray of 3 Compartments
- 1 x Drape blue wrap 40 x 40cm sterile field
- 1 x Drape blue napkin 40 x 40cm
- 4 x Gauze Swab 3 x 3 8ply
- 6 x Absorbent Cotton Balls
- 2 x Glove S



Biopsy Dressing Set M

- 1 x Plastic Tray
- 2 x Plastic Cup 30ml
- 1 x Fenestrated Drape 60 x 60cm
- 1 x Drape blue wrap 75 x 75cm
- 4 x Gauze Swab 3 x 3 8ply
- 2 x Sponge handle
- 1 x Plastic forceps
- 1 x Blade #11



Biopsy Dressing Set S

- 1 x Plastic Tray
- 2 x Plastic Cup 30ml
- 1 x Fenestrated Drape 60 x 60cm
- 1 x Drape blue wrap 75 x 75cm
- 4 x Gauze Swab 3 x 3 8ply
- 6 x Absorbent Cotton Balls
- 2 x Plastic forceps
- 1 x Blade #11

