



ENOUSAFE



# RADIO ANGIOGRAPHY PACK

Disposable Radio Angiography Pack is optimise efficiency for cath lab cases, Body & neuro - interventional cases and or cases.

This pack is sterile set which help to reduce cost, setup time and safe for both of physician and patient.

The component in pack include radial & femoral angiography drape,surgical gown,table cover,wrapping drape,gauze,hand towel,sponge handle,plastic tray&plastic cup,band bag,screen cover, remote control,blade and color sticker.



DISPOSABLE RADIO - ANGIOGRAPHY PACK

# RADIO - ANGIOGRAPHY PACK

## 1. Radial & Femoral Angiography Drape

Radial and Femoral fenestrated drape made from Trilaminare material.

- Dual Clear Side Panels 50cm width
- Radial & Femoral fenestration Ø12cm
- SMS full Absorbent Reinforcement 140cm x 370cm



## 2. Surgical gown



Gown SMS made from Polypropylene non-woven textile of 3 layers Spunbond Meltblown Spunbond (SMS)

- Reinforce in chest and sleeves excellent physical barrier to contaminant and fluids.
- The Polypropylene fibres have advantage of building barrier against bacteria and repelling blood and other liquids. With good tensile strength, high filtration, lint-free, breathable.

## 3. Table cover , Wrapping Drape

- Made from dual-layer material with its absorbing and liquid trapping property, quickly absorbs any type of liquid that comes onto the cloth and PE film creating the bottom layer prevents absolute liquid penetration. Its high barrier against bacteria, impermeability to blood and other liquids, lint-free, good drapability.



- Made of natural fiber-based materials, each sheet is biodegradable and no emit odors, toxic waste or other residues in an incinerator.

## 4. Hand Towel & Gauze

- Spunlace non-woven consist of fleece material, structure due to its skin comfort with its soft textile, smooth surface lint free, high wearing comfort, breathable, tear-proof, against bacteria.



- Gauze 10x10cm 12ply, made from pure 100% cotton, high absorbent, softness

## 5. Plastic Tray & Plastic Cups

- Plastic Tray made from food grade PVC.
- Plastic Cups and kidney dish made from Polypropylene (PP) material.



## 6. Band Bag , Screen Cover , Remote Control


- Band Bag latex-free, elastic at opening, 3s' thickness



- Screen Cover light cover square 100x120 cm, 3s' thickness



- Transparent PE Plastic for Cover Medical Protective Equipment) , 3s' thickness + Rubber ring



## 7. Sponge handle , Blade , Color Sticker

- Plastic handle medical sponge length 16cm



- Disposable plastic handle blade No.11



- Water proofs (White,Green,Yellow, Blue,Red)



Component in pack include :

- |  |                                   |
|--|-----------------------------------|
| 1 (EA) Radial & Femoral Angiography Drape size 240cm x 370cm | 2 (EA) Plastic Cup 50cc           |
| 2 (EA) Surgical gown size L                                  | 2 (EA) Plastic Cup 250cc          |
| 1 (EA) Surgical gown size XL                                 | 1 (EA) Plastic Cup 1000cc         |
| 1 (EA) Table cover 140cm x 230cm                             | 1 (EA) Kidney dish 700cc          |
| 1 (EA) Wrapping Drape 120cm x 120cm                          | 3 (EA) Band Bag Ø70cm             |
| 2 (EA) Hand Towel 30cm x 40cm                                | 1 (EA) Screen Cover 100cm x 120cm |
| 5 (EA) Hand Towel 55cm x 78cm                                | 1 (EA) Remote Control 10cm x 26cm |
| 30 (EA) Gauze 10cm x 10cm 12ply unfold                       | 4 (EA) Sponge handle 16cm         |
| 1 (EA) Plastic Tray 14 inch x 14 inch                        | 1 (EA) Blade #11                  |
|  | 1 (EA) Color Sticker              |

